


# Luwa® JM Sterile Air Filter



- Efficiency > 99.99995 %
- Compact and robust construction
- Extreme security due to two stage filter element
- Low pressure drop
- Compact stainless steel case, resistant to external pressure of up to 10 bar
- Individually tested and leak free

**Luwa**

# Luwa JM Sterile Air Filter

Luwa JM sterile air filter is a cylindrical HEPA filter in compact metal case fitted with 1/2 inch gas male connections at each end. It is designed to filter small air flows at pressure up to 10 bar and it offers, due to its two stage filter element, an extremely high efficiency and security against leakage.

## Fields of application

Filtration of compressed air and gases in:

- Medical technology
- Food and beverage industry
- Pharmaceutical industry
- Clean room technology

The Luwa JM sterile air filter is used for:

- Anaesthetic and respiratory apparatus
- Ventilation of autoclaves after release of vacuum
- Particle counters, as zero indicator filter
- Filtration of pressured air for regulating and control units
- Cleaning and drying of surfaces by means of compressed air
- Storage tanks and reactors as breather filter
- Others

## Construction and materials

### Casing:

Electropolished S 304 stainless steel (1.4301)

### Filter medium:

Two identical panels of pleated micro-fine glass paper in series. Each panel being individually sealed and supported and having protection on each side from mechanical damage.

### Sealant:


Polyurethane

### End connections:

S 304 stainless steel 1/2 inch gas male connectors welded to the casing and fitted with spanner flats. Factory fitted with plastic caps.

### Sterilisation

The JM filter is not sterilised in the factory. A 20 minute autoclave 134 °C would ensure complete sterilisation of the filter. The filter would require subsequent drying at less than 90 °C until it is within 2 grammes of its original weight. Sterilization might also be achieved by spray disinfection. Weight control before



## Technical data

Type	JM
Part number	424 10 0101
Efficiency as per EN 1822 (MPPS <sup>1</sup> )-DEHS aerosol)	> 99.99995 %
Rated air flow	10.0 m <sup>3</sup> /h
Initial pressure drop ( $\Delta P_A$ )	$\leq 1000$ Pa
Recommended final pressure drop	$\Delta P_A + 600$ Pa
Maximum final pressure drop	10000 Pa/0.1 bar
Maximum operating pressure <sup>2)</sup> (across casing)	10 bar
Maximum continuous operating temperature	110 °C
Maximum sterilization temperature <sup>3)</sup>	140 °C
Maximum relative humidity	< 100 %
Approximate weight	0.6 kg

usage (see above) is then also required. The validation of the filters function after sterilization has to be performed by its user.

### Installation

The JM filter is fitted with 1/2" gas threaded male stubs complete with spanner flats which must be used in tightening the filters in a pipe line system to prevent damage of the housing. The threads should be sealed with PTFE or similar tape or by means of compression rubber seals. The air flow direction arrow on the housing must be adhered to. For installation in pipe networks with male threading use ordinary joints (e.g. from Serto, +GF+, etc.).

### Tests

Each filter is tested to DIN 24184 (part 4.2) to ensure absence of leaks. The casing tightness is individually tested to 6 bar. A tightness test to 16 bar is avail-

able on request. The burst pressure of the unit has been verified on type tests to be greater than 110 bar.

### Remarks

- <sup>1)</sup> MPPS (**M**ost **P**enetrating **P**article **S**ize): the particle diameter at which the separation efficiency reaches its minimum – usually in the range 0.1–0.2  $\mu\text{m}$
- <sup>2)</sup> Temperatures of 20 °C
- <sup>3)</sup> For short time in saturated vapour (see "sterilisation")

# Luwa JM Sterile Air Filter



## Technical data

**Initial pressure drop ( $\Delta P_F$ )**  
as a function of air flow ( $\dot{V}_L$ ).

Valid for ambient conditions:

$t_L = 20\text{ °C}$   
 $\rho_L = 1.15\text{ kg/m}^3$   
 $p_{atm} = 0.966\text{ bar}$

for connection pipe with  $\varnothing 14\text{ mm}$  inner diameter.



### Initial pressure drop ( $\Delta P_F$ )

as a function of air flow ( $\dot{V}_N$ ) and absolute pressure ( $P_1$ ) at  $T_1 = 273\text{ K}$ .

Calculation of  $\dot{V}_N$ :

$$\dot{V}_N = \dot{V}_1 \cdot \frac{\rho_1}{\rho_N} [\text{m}^3/\text{h}]$$

$$\rho_1 = \rho_{tr} (1 - 0.377 \cdot \varphi \cdot \frac{P_p}{P_1}) [\text{kg/m}^3]$$

$$\rho_{tr} = 387.75 \cdot \frac{P_1}{T_1} [\text{kg/m}^3]$$

with:

$t_N = 0\text{ °C}$   
 $P_N = 1.013\text{ bar}$   
 $\rho_N = 1.293\text{ kg/m}^3$  } Normal conditions

$\dot{V}_1$  = air flow as a function of  $(P_1, T_1, \rho_1)$

$\rho_1$  = air density at condition ① [ $\text{kg/m}^3$ ]

$\rho_{tr}$  = dry air density at condition ① [ $\text{kg/m}^3$ ]

$\varphi$  = relative humidity

$P_p$  = partial vapour pressure [bar]

$P_1$  = absolute pressure [bar]

$T_1$  = absolute temperature [K]

